


# How to use thermometer

1. Remove the 98.6°F sticker from the thermometer's display.
2. Find the thermometer's button at the top and click it once to turn it on.
3. Wait for APP to appear on the display
4. With "APP" showing on the display, VERY QUICKLY double-click the same button. (Think the same speed as double-clicking a mouse or retractable pen.)

When the double-click is successful, you'll hear a short high-pitched beep and the screen will change from "APP" to three lines "- - -". This means the thermometer is now ready to take a reading. The thermometer can now be used like any other digital thermometer from here on, and you will not need to double-click the button again in the future.

5. to turn off thermometer press the button and hold down.



# Como usar su termómetro

1. Remueva el sticker de “98.6 F” de la pantalla del termómetro
2. Encuentre el botón hacia arriba en su termómetro y presionarlo una vez para prenderlo
3. Espere que aparezca “APP” en la pantalla
4. Con “APP” en la pantalla, RÁPIDO presione el mismo botón dos veces  
(Igual como presiona un lapicero retractable )

Si presiono dos veces exitosamente, escuchará un sonido corto y la pantalla cambiará de “APP” a tres líneas “- - -.” Esto significa que su termómetro esta listo para tomar la temperatura. El termómetro de aquí en adelante esta listo para usar y no necesita presionar dos veces en el futuro.

5. Para apagar su termómetro presione el botón por unos segundos

